

LockeT

Intuitive Suture Retention Device

The LockeT device offers a cost-effective solution for site closure with suture retention, ensuring efficient hemostasis while prioritizing patient comfort. Compact and user friendly, LockeT simplifies the closure process, making it easy to monitor, adjust and remove as needed.



Simple design,
easy to use



Patient
Comfort



Cost Effective



Suitable for a wide
range of catheter
sizes up to 27F




Quick
Hemostasis




Transparent
design for easy
monitoring

LockeT


Intuitive Suture Retention Device



Step 1: Create a standard figure of eight and remove sheaths.



Step 2: Tighten sutures, cleanse site and place sutures within LockeT.



Step 3: Rotate the handle 90 degrees.

ORDERING INFORMATION

MODEL NUMBER	PRODUCT DESCRIPTION	QUANTITY
9300	Suture Retention Device	1 Box (10 devices/box)

LockeT is a Class 1 sterile product registered with the FDA under Device ID D496327

Refer to IFU for complete Indications for Use, contraindications, warnings, and precautions.